


**A12Q MASSAGE CHAIR**  
**Operation instructions**

SAVE THESE INSTRUCTIONS



Thanks for purchasing this product. Please read this manual instruction carefully in order to know how to use this product correctly before operation. Please pay more attention to some important safety warning and keep this manual book properly for future use.

Note: Our Company reserves the right to revise the design and description of this product without any further notice. Please refer to the actual color of this product.

## CONTENTS

---

|                                 |    |
|---------------------------------|----|
| Safety and maintenance          | 02 |
| Name and function of components | 04 |
| Function & Efficacy             | 05 |
| Method of usage                 | 06 |
| Product specification           | 17 |

# IMPORTANT SAFETY WARNING.

---

- Don't allow children to touch moveable part of this product.
- Please use well-grounding power supply suited to this product.
- Please pull out the plug after use or before cleaning to avoid injury or damage to the product.
- Please operate this product according to this manual instruction.
- Don't use accessories which are not recommended.
- Don't use this product outdoors.
- Please read this manual instruction carefully before operation.
- Any other usage that is not listed in this manual instruction is forbidden.
- 15 minutes of usage each time is recommended.
- Please don't use it if the leather is broken or damaged.
- Please don't use it if the cover or leather is damaged.
- Please don't use it if the radiator hole is covered or blocked.
- Please don't drop anything into this product.
- Please don't fall asleep while using this product.
- Please don't use it if drunk or feel unwell.
- Please don't use it within one hour after having meal.
- Please don't make the massage function too strong for avoiding injury.

## 2. Environment for usage.

- Please don't use it under high temperature and moist environment such as bathroom.
- Please don't use it immediately while the environmental temperature changes sharply.
- Please don't use it under heavy-dusty or caustic environment
- Please don't use it where there is not enough space or not good ventilation.

## 3. People who are not suitable to use this product.

- People who are suffering from osteoporosis are not suitable to use this product.
- People who are suffering from heart disease or carrying electronic medical equipment such as heart pacemaker are not suitable to use this product.
- People who are suffering from fever are not suitable to use this product.

# Safety and maintenance

---

- Pregnant women or women that are in menses period are not suitable to use this product.
- People who got injured or skin disease are not suitable to use this product.
- Children under twelve or people who are not able to operate this product are not suitable to use this product.
- People who are told by doctor to have rest or people who feel unwell are not suitable to use this product.
- People whose bodies are wet are forbidden to use this product.
- People whose is sensitive to heat

## 4. Safety.

- Check the voltage if it is suitable to the specification of this product.
- Don't use wet hand to pull the plug.
- Don't make water spatter to this product to avoid electric shock or cause damage to this product.
- Don't damage the wires or change the circuit of this product.
- Don't use wet cloth to clean the electric parts such as switch and plug.
- Be away from this product under power cut state to avoid injury if the power comeback suddenly.
- Stop using this product while it works abnormal, and please consult local agent immediately.
- Stop using this product if you feel unwell, and please consult health care doctor.


## 5. Maintenance.

- This product should be maintained by certain appointed agents, users are prohibited to disassemble or maintain by themselves.
- Please don't forget to cut the power after usage.
- Don't use this product if the socket looses
- If this product will be left unused for long, please curl up the wires and reserve this product in dry and dustless environment.
- Don't reserve this product under high temperature or near fire and avoid longtime direct sunshine.
- Please clean this product by dry cloth. Don't use thinner, benzene or alcohol.
- The mechanical components of this product are specially designed and made, it is not necessary to be specially maintained.
- Don't use edged res to thrust this product.
- Don't roll or pull this product above uneven ground, it should be lifted before moving.
- Please use it intermittently; don't make this product work continuously for long.

## 6. Solution of ordinary malfunctions.

- It is nature movement voice if the motor sounds while in usage.
- If the controller can't work normally, please check the plug and socket if they are connected strongly and make sure the switch is turned on.
- If the rated working time is over, the product will automatically activate the power-off button; If this product works continuously for long, the temperature protector will make the product automatically power-off, it should be used again after half an hour's rest.

## Name and function of components



## Function & Efficacy

### Function

- Massage mechanism is four-wheel driven, muted design can move up and down.
- Shiatsu, kneading, flapping, knocking, massage and rolling six emulational massage models.
- Automatic check the upper body massage points, shoulder position check and adjustment.
- Five characteristic massage function (comfort, body care, ease, activate Relax).
- Upper body automatic massage (shoulder, back, waist overall, partial and Fixed six choices).
- Upper body manual massage (six choices: shoulder, back, waist overall, partial, fixed), it is designed with five speed levels for each massage choice. under the state of flapping, shiatsu and Knocking, the width of the massage mechanism can be changed in three Degrees: wide, narrow and medium.
- Under fixed state, the massage mechanism can move up and down to reach accurate position.
- Lower body airbag massage function: equipped with 33 airbags, three modes and three levels of speed. (Strong, medium and weak).
- Arms airbag massage function: equipped with 36 airbags.
- Vibration massage function: set in the back and buttock position, there are three choice: strong, weak and off.
- Music display and volume adjustment function; previous song, next song choice, r.
- Music accompaniment function: back massage speed and vibration changes with the high frequency, buttock vibration strength changes with the middle frequency.
- Under the music accompaniment state, the controller in the wire can be used to adjust the massage strength.
- Automatic recline and decline function: it is designed with modes that the calf rest and backrest recline and decline separately, lie down.
- Recover of the massage mechanism when turning off the chair.
- It is equipped with one heater which use five Xiu jade as the infrared emission source, the temperature can be adjusted between 40°C ~ 70°C.
- It is equipped with seven color LCD, MP3 displayer, earphone and mini tea Table.

### Efficacy

Dispel fatigue.  
Accelerate blood circulation.  
Alleviate muscle fatigue.  
Relax tense muscle.  
Alleviate nerve aches.  
Alleviate muscle aches.

## Method of usage

### 1.The usage of the control

- on/off massage

Press power button to start massage. During the massage period, press the power button and it will stop, the calf rest will return to its original position( move down to the bottom).


- Neck position adjustment: press up and down button.

### 2.Characteristic automatic massage


Comfort: massage whole upper body and meantime start all other functions (Lower body airbag and vibration massage) to relax muscle. Set the working Time to be 8 minutes. Mode display window shows (01).



Body care: besides the whole upper body massage function, it can massage respectively shoulder, back and waist. Meantime start all other functions to accelerate blood circulation. Set the working time to be 15 minutes. Mode display window shows (02).




Ease: massage shoulder, back and waist, meantime start all other function to relax aches in the shoulder, back and waist. Set the working time to be 30 minutes. Mode display window shows (03).




Activate: deeply massage whole upper body, shoulder, back and waist, meantime start all other functions to relax whole body. Set the working

## Method of usage

Time to be 38 minutes. Mode display window shows (05).



Relax: total 16 massage functions. Divide the upper body to be 16 massage points, do deep massage according to your choice to make the massage more accurate. Relax the aches of the muscle in the massage points and stimulate points beside the acantha.




Awake: automatically start massage function. During the massage period, press (awake) button then all massage function will stop and the timing display window shows (01), then press again the (awake) button to set the automatic start time. The automatic start time can also be set when the massage chair is turned off. The max timing is (08) (that is 8 hours). After setting the automatic start time, you may lie on the chair to sleep or relax and you won't be worry about the time. Because half an hour before the automatic start time ends, the chair will be start automatically to remind you that time is up.

Timing display window shows the automatic start time.

Mode display window shows (99)


E.g.: press (awake) button continually three times, then the timing display window shows (03) which means the automatic start time is 3 hours, and the mode display window shows (99) which means the awake function is turn on. From the moment you finished setting the time, two hours and a half later the massage function will start automatically.




## Method of usage

3.Auto---timing is 15 minutes.


- Whole upper body automatism---shows (F0).




- Partial automatism---shows (F1). Press up and down button to move the mechanism to the required position.



- Fixed automatism---shows (F2) Press up and down button to move the mechanism to the required position.




- Shoulder automatism---shows (F3).




## Method of usage

- Back automatism---shows (F4).



- Waist automatism---shows (F5).



4.Manual function.

- Timing: it can only be set under manual massage function; it is designed with three levels: 15 minutes, 20 minutes and 30 minutes.


- Manual: total 6 manual choices: fixed, partial, overall, shoulder, back and waist shows respectively P0, P1, P2, P3, P4 and P5.

- Massage methods:


Flapping, knocking, massage, rolling, shiatsu and kneading.



- Width: adjustment of the width of massage mechanism. It is designed with three levels: wide, middle and narrow



Wide



Middle



Narrow

- Speed: adjustment of the massage speed. The speed can't be changed under the shiatsu mode, and under other modes, it can be changed with five levels.



## Method of usage

- Back vibration: back vibration massage strength and on/off button.



- Buttock vibration: buttock vibration strength and on/off button.



- Arms: arms airbag strength and on/off button.



- Lower body: lower body airbag on/off button.



- Strength: lower body airbag strength button, three levels: strong, medium and weak.



- Mode: lower body airbag modes choice. There are three modes: Lower body, Foot & Calf, Buttocks.



- Handed heat therapy: on/off button of the handed heat therapy. Default temperature is 40°C. (there are seven levels available: 40°C、45°C、50°C、55°C、60°C、65°C、70°C).




- Handed heat therapy(+): increase the temperature of the handed heater. It will be increased by one level upon each press, the highest temperature is 70°C.



- Handed heat therapy(-): decrease the temperature of the handed heater. It will be decreased by one level upon each press. The lowest temperature is 40°C.

## Method of usage



- Background color: the background color changes automatically with the massage methods.

- Color adjustment: set manually the background color. There are seven colors: red, green, blue, purple, orange, white and cyan.


### 5.MP3 display

- Function: turn on, turn off, previous song, next song, synchronization.



- Music accompaniment function: (1. Back massage speed and vibration strength changes with the high frequency. 2. Buttock vibration strength changes with the medium frequency.)

- 6.Calf rest and backrest recline and decline function. (It is designed with calf rest recline and decline, backrest recline and decline, sit and lie down).



## Method of usage

---

### 7.How to download mp3 music :


1. Turn off the power of massage chair .
2. Pull out the u-disk from the MP3 player .
3. Insert the u-disk into USB port of PC (Windows Xp、 Windows 2000、 WindowsNT.)
4. Download your favorite MP3 format music into U-disk .if fails , please firstly format the disk.
5. Delete the USB Mass Storage Device on the PC after complete download .
6. Pull out the U-disk from PC
7. Insert U-disk to the USB port of MP3 player
8. Now MP3 music was renewed . you may turn on the power of massage chair .

### Note :

1. If fails to play MP3 music ,please firstly format the U-disk on PC before download.Then follow the instruction to insert the U-disk into USB port of MP3 player.
2. Except the U-disk ,do not remove or loose the MP3 player .

## Method of usage

---



( 1 ) MP3 player

( 2 ) U-disk front cover

( 3 ) U-disk


## Method of usage

### 8. THE PADDING AND BACK CUSHION ASSEMBLING AND OPERATION

Using the cushion which can reduce the neck & shoulder kneading massage of the pillow, can determine according to your demands whether to use the pillow cushion or not (recommend using the pillow cushion).


backs of cushions connect back through zipper (1) pillow cushion and back through its spike joining to glue cushion (2)

**Warn:** Can't raise the back cushion to use, in case that cause the injury to you.




### 9. INSTALLATION AND USE OF THE CONTROLLER AND TEA TABLE TRAY

Install the controller shelf on the massage chair (1) Insert the controller shelf in the controller shelf (2) The control can rotate freely and vertically (3) The control can rotate freely and horizontally (4) Install the tea table tray on the massage chair (5) Put the tea table tray on the shelf (6)



## Method of usage




### 10. PROHIBITED PROCEEDING:

- (1) When the arms airbag groupware connect with the armrest, the spindle can't be added with lube.
- (2) It is forbidden to stab the airbag by needle or other sharp things.
- (3) It is forbidden to wobble the airbag groupware in case to destroy relevant components.

### Usage:


- (1) When the airbag groupware spindle is inserted into the hole half length and the chair is turned on, the induction equipment in the hole receive signal and will start feeding the air; when the groupware spindle is inserted into the hole completely, there will be no air coming out of the joint and the airbag will start working.
- (2) When the airbag groupware is taken off but not completely, it is common that there will be air coming out of the joint.
- (3) The hand-wheel will be loose if turned to left side, and be tight if turned to right side; it can be taken off by turning to the left side continuously.
- (4) When take off the groupware, please check if it moves forward and backward freely.
- (5) In order to avoid colliding with the airbag groupware, the controller should be put forward by 20° as shown in picture (1) and (2).



## Method of usage

(6) the process of installation of arm massage


- 1) Get down the rubber stopple.
- 2) Put the rubber stopple into the accessory case.
- 3) Adjust the both bars to the corresponding holes, and push the machine hard, rotate the handheld clockwise and tighten it up.



### 11. FLOOR PROTECTION

Put the heavy massage chairs on the wood floor for a long time may damage the floor, so please place a carpet or other object to avoid this damage.

Caution: Do not move the chairs with its wheels on wood floor, rough ground or in a narrow space. Two people are required to lift the chair.



## Product specification

Model: A12Q

Description: Massage chair

Voltage: 110-120V~ 60Hz ☐  
220-240V~ 50Hz/60Hz ☐

Power: 220W

Safety design: Class I

Rated working time: 15 minutes